


NeuroShield

FOR NERVE REPAIR



TECHNICAL GUIDE




DEFINITION

NeuroShield™ is a chitosan membrane to provide a non-constricting protection for peripheral nerves.

NeuroShield™ is designed to be an interface between the nerve and the surrounding tissue for uses to treat nerve injuries. When hydrated, NeuroShield™ is easy to handle, soft, pliable, nonfriable, porous and transparent.

NeuroShield™ is provided as a sterile, non-pyrogenic rectangular sheet in the size of 40 x 30 mm and is intended for single use.

NeuroShield™ can easily be placed over the injured nerve, and can be sutured if necessary. Furthermore, the device can be trimmed or shaped to the appropriate size to fit the nerve tissue to be treated.



- Porous membrane
- Biodegradable
- Flexible and pliable
- Clinically proven



ORDERING INFORMATION

NeuroShield™ (40 x 30 mm) Code NS4030



Monarch Bioimplants GmbH
Platz 4
CH-6039 Root D4
Switzerland

+41 (0)41 4552261 / +41 (0)41 4552262

Fax: +41 (0)41 4552264

mail@monarchbioimplants.com

www.monarchbioimplants.com

