




M·CATH
FLEXY

Track and Push

M-CATH Flexy


Microcatheter

- Reinforced tip
- Lengths 135 cm and 150 cm
- Tapered balanced braided shaft
- Durable hydrophilic coating
- Enhanced torsional capacity



M-CATH Flexy

Microcatheter




Reinforced tip



Great flexibility



Comfortable handling



Ordering Information/ Technical Information

Reference	Usable length	Tip entry profile	Dist. Outer Diameter	Prox. Outer Diameter	Guiding Catheter Compatibility	Max Guidewire Size
820135732	135 cm	1.65 F 0.55 mm	2.46 F 0.82 mm	2.79 F 0.93 mm	5 F 1.67 mm	0.014" 0.36 mm
820150732	150 cm	1.65 F 0.55 mm	2.46 F 0.82 mm	2.79 F 0.93 mm	5 F 1.67 mm	0.014" 0.36 mm

Track and Push